

KALSTEIN® Phototherapy Lamp

Model YR02194



Characteristics:

- ✓ LED fluorescent lamp as the phototherapy radiating elements
- ✓ Angle of lamp head can be adjusted horizontally.
- ✓ Height of lamp head can be adjusted.
- ✓ Stainless steel support stand can be adjusted.
- ✓ Wheels can be locked
- ✓ Timer to record the phototherapy operation.
- ✓ Use in incubator, irradiant heat bassinets, and any kind of bassinet.

Application

The YR02194 unit like YR02195 unit with blue fluorescent tube has independent framework can use with incubator, irradiant heat bassinets, and any kind of bassinet in family care. Widely applicable in all hospitals and clinics and optimize for jaundiced baby.

Phototherapy Lamp, LED panel, mobile w/ braking castors.

The newborn patients with hyperbilirubinemia all need best care for threatening jaundice. Our line of phototherapy products supply you the reliability and stability intensive phototherapy by LED light for your baby patient. The unit has a long-life that provides 20,000 to 0,000 hours of light therapy.

General Features:

- ✓ Big area lighting to ensure sufficient treatment effect with uniform intensity
- ✓ Fan-free design to eliminate noise and disturbance to laminar flow
- ✓ Innovative LED (NEW) and tube fluorescent lamp (conventional) available
- ✓ Models with fixed or adjustable working angle are available
- ✓ Installed on stretchable stand-pole or integrated over incubator with 4-leg sucker device
- ✓ 99 levels of adjustable light intensity for individual phototherapy need
- ✓ Time countdown setting to avoid overtreatment
- ✓ Ergonomic design for easy access to baby cradle, infant warmer or incubator.

Model	YR02194
Function	
Irradiation wavelength range	400-550nm
Adjustable Irradiation intensity (max)	3000 uw/cm2
Wavelength of wave crest	460nm
Adjustable Irradiance from 1 to 99 levels	√
Timing	√
Fixed time	√
Time totalizer	√
Display range (hour)	0-9999h
Structural parts	
Vertical Height Adjustable	1100-1580 mm
Useful service life	5000h
Nominal blue light life	20000h
Dimension	500x600x (1355-1705)mm
Electric parts	
Power supply	AC(220 ±22)V (50±1) Hz
Power consumption	≤50W
Ambient temperature	10-40
Ambient relative humidity	≤80%
Motor	AC



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663210023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S

