


SCBI for STEAM, Plasma & Ethylene-Oxide



SAP / Monitoring

Quickly monitoring with SCBI

The monitoring of sterilization condition at preferred places with the help of Self Contained Biological Indicators (or SCBI) is no longer practice of every day in most European Hospitals. Whereas outside Europe the SCBI is widely used for daily routine checks for Steam sterilization in Hospitals.

REF: 107.001.0100	SCBI for STEAM, pop. 1x10E6	Unit: 100 pieces
REF: 107.002.0100	SCBI for Plasma, pop. 1x10E6	Unit: 100 pieces
REF: 107.003.0100	SCBI for EO, pop. 1x10E6	Unit: 100 pieces
REF: 120.001.0001	Dual Incubator for SCBI	Unit: 001 piece

Characteristics:

- ◆ Biological Indicators conform to ISO 11138
- ◆ SCBI EO: Bacillus Atrophaeus ATCC 9372 1x10E6 spores per vial
- ◆ SCBI STEAM/Plasma: Geobacillus Stearotherophilus ATCC 7953 1x10E6 spores per vial
- ◆ Color change: Green to yellow (EO) and purple to yellow (STEAM)
- ◆ Digital Incubator (220V) for 37°C / 48 hrs (EO) and 56-59°C / 24 hrs (STEAM & Plasma)

